


EXPERIMENT NO. 4

DETERMINE THE DISPERSIVE POWER OF Na D-LINE USING GRATING.

INTRODUCTION: Departure from the expectations of ray optics are broadly classed under the term diffraction. A diffracted image thus means an image off the line of sight due to deviation of light from the straight line path. A series of such images at increasing angles of deviation are called as 1st, 2nd, 3rd order images respectively shown as P_1, P_2, P_3 etc in Fig. (1). P_0 is the direct image (central order) where all diffracted rays are meeting in phase.



The condition to obtain P_1, P_2 etc. images (i.e. points of maximum intensity) is that the path difference at a given point on the screen, between the two rays coming from the corresponding points of the slits is an integral multiple of, i.e. the wavelength of incident light. Clearly from the Fig. (1), this condition will be represented by the equation.

$$(a + b) \sin \theta = n \lambda$$

$$\text{Dispersive power} = d\theta / d\lambda = n / (a + b) \cos \theta$$

where $(a + b)$ is grating element

a = width of a transparency

b = width of an opacity

n = order of the spectrum

θ = angle of diffraction

Measure θ as explained below and make use of the above equation to find the dispersive power of Na D-line.

EXPERIMENTAL DETAILS:

- (i) **SETTING OF SPECTROMETER FOR PARALLEL RAYS:** Point the telescope towards a distant object with a marked outline e.g. an edge of a wall. Seeing through the telescope first make the cross-wire distinctly clear adjusting the eye-piece tube.